COMPRESSANA

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Instructions for putting on support stockings

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
Instructions for putting on

en


1 Turn sock inside out as far down as the heel. The foot part should remain inside the sock. Stretch the opening for the foot apart with both thumbs/hands and pull the sock over the foot up to the heel.

2 a) Then turn over the sock section which has been turned inside out by pulling it up or pushing in the leg. This will produce a double crease in the middle of the foot.

b) Reach into the double crease with two fingers and lift the sock over the heel.

3 Hold the sock with both hands and place it on the lower leg by turning it up in stages without snagging. Smooth any creases which are formed by pulling upwards.

Ideally put the socks on after you get up in the morning. Make sure that there are no creases, that the rows of stitches are straight and that the toes are not constricted. Do not turn the socks fully inside out in order to prevent any constriction. Wearing regularly is recommended.

COMPRESSANA GRIP rubber gloves make it much easier to put on the socks. A special surface structure ensures that the knitted fabric can be pulled upwards in stages. It is just as easy to correct the fit of the socks later on and remove any creases and twists. In addition, when you put the socks on, this prevents any damage (e.g. from long finger nails or jewellery with sharp edges).

Thanks to the smooth and slidable material of the **COMPRESSANA EASY** pull-on aid, the sock can be pulled onto the leg without any friction. EASY is suitable for both open and closed toes.

COMPRESSANA products can be used in combination with **COMPRESSANA SKIN CARE preparations**. After wearing these socks, make sure you perform regenerative skin care in order to prevent any dryness and itching.

The **COMPRESSANA vein ointment** cares for and protects your legs and is fully absorbed into the skin. It is therefore possible to use it before and after wearing socks.

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Compressana GmbH • Böhmerwaldstr. 3