TAPE CORRECTION SOX

TYP _____
**HALLUX
VALGUS**

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Japanese tabi socks with separate toe areas or special socks with integrated non-adhesive and skin-friendly tapes

- + Combat misalignment and correct malpositioning of toes or foot, depending on the model.
- + Stabilise the toes, mobilise joints, promote a natural rolling movement of the foot, reduce pain and can prevent arthrosis.
- + Suitable for wearing all day and gentle on the skin like normal socks/stockings; discreet can be worn without restrictions and suitable for the workplace, leisure and when exercising
- + mobilisation = relief + regeneration

Tapes are integrated into the inside of the Taping Socks. They vary according to the model and aim to counteract typical malpositioning.

The Taping Socks combine the benefits of conservative methods of treatment, such as physiotherapy, with kinesiotaping to counteract malpositioning of the foot or toes. They are free from side-effects. They can be used long term for problems that have functional causes. These include: bunions, flat feet, splayfeet, hammer toes.

Application areas

Preventative: The COMPRESSANA TAPE SOX can be worn as a precautionary measure if you have an existing or developing misalignment. It is possible and generally unproblematic to combine them with inserts.

Post-operatively: After an operation for bunions or hammer toes, taping socks can help maintain a successful operation result by countering local misalignment.

The TAPE SOX can also be worn at night for gentle, passive correction. The corrective socks offer relief and support for flat feet and splay feet.

TYPE
HALLUX
VALGUS

Bunions (Hallux valgus)


Bunions are the most commonly occurring and treated cause of toe misalignment. Bunions cause the big toe to be malpositioned towards the outer edge of the toe. Women are more prone to this type of foot deformation than men. If left untreated, bunions can cause hammer toes.

TYPE
HAMMER
TOE

Hammer toes / claw toes


A hammer toe is a toe (usually the second toe) that is stiffened and bent upwards, resembling the shape of a hammer. This malpositioning becomes increasingly painful. Again, women are about five times more likely to suffer hammer toes than men. In the early stages the toes can often be moved back to their normal position. As it progresses, the toes become stiffened and arthrosis develops. Hammer toes and bunions often occur together.

TYPE
FALLEN
ARCH

Flat feet, splay feet (Fallen arch)

If the heel turns inwards and forwards when weight is placed on it, the longitudinal arch becomes flattened. This is termed low arch and flat foot. This causes instability of the outer edge of the heel, making it turn inward. When the body weight is on the foot, there is also pressure on the longitudinal arch. If left untreated, this can cause flat foot. If the shoe sole is worn down towards the inside in the heel and forefoot areas, this indicates overpronation due to fallen arches/ flat foot.

Use

The best results are achieved from COMPRESSANA TAPE SOX if worn in shoes that are not too tight, with flat heels (< 4 cm) and enough room for the toes. The shoes must not be too small, as this pushes the toes out of the correct position, which can cause bunions, as well as hammer/claw toes. It is important that you slowly increase the length of time that you wear the corrective socks, or wear them first overnight, so that the toe joints can gradually get used to the new mobility. In the course of a week you can increase the length of time you wear COMPRESSANA TAPE SOX until it is for the whole day.

Instructions for use

Bunions / hammer toes: When putting on the socks, please note the separate toe areas. First, gather up the sock and position the big toe in the correct compartment. Do the same for hammer toe corrective socks, putting the second toe in the little box.

The taping band for the hammer toe must be pulled over the middle bone towards the base joint, so that the pull band is over the raised/bent part of the toe (see image). The hammer toe tape must not be in front of the raised part as there is a significant risk that this could worsen the misalignment.

Next, position the other toes in the remaining toe boxes and then pull the Taping Sock towards the heel, making sure the pull ties are in the right place. To ensure correct functioning, the taping bands must be under tension with the toes fully in the toe boxes and as taut and crease-free as possible. If the sock is too loose, the corrective function is limited. In this case please choose a smaller size.

Pronation Control: When putting the socks on, make sure you have the correct socks for right and left. The words "Taping Socks" should be on the outer side of the foot.

The taping bands are divided into two function areas. The heel band goes on the outside, from the base of the Achilles tendon, under the outer ankle to the foot sole. The upper edge of the taping band lightly touches the outer ankle ①.

The front taping band goes from the foot sole at a 45 degree angle to the ball of the foot at the back of the foot. The front edge lies behind the little toes ②.

On the inner side the band runs from the start of the Achilles tendon under the inner ankle towards the foot sole. Here the upper edge of the taping band lightly touches the inner ankle ③. Pulling the taping band stabilises the heel in the correct position and raises the inner longitudinal arch ④.

The front taping band on the inner side stabilises the first metatarsal bone towards the ground ⑤. Please use the two finger technique to check the position. Check that the lower leg is at a 90° angle to the foot by putting your index and middle fingers in the instep ⑥. If the back edge of the taping band is below your middle finger, it is correctly positioned ⑦.

Wearing time = treatment time: You can wear taping socks inside slippers, as well as outdoors and sports shoes. They can be worn for long times. As taping socks they are suitable for wearing with sandals or without shoes.

The familiarisation period

It's important that during the familiarisation period you also strengthen your muscle structure. There are various exercises that help to relax and strengthen your muscles.

An example for bunions:

Exercise to strengthen the big toe abductor

- Put your foot on the ground.
- Put more weight on the outer edge, particularly the outer ball than the inner side
- Brace the big toe inwards against the floor
- Hold the tension for 10 seconds
- Relax the foot
- Repeat the exercise 3–4 times

Try to do the movement in a controlled way, consciously using the muscles. Often the nerve impulse (conscious action) has been lost. So first it needs to be controlled, then strengthened as above.

Your physiotherapist will give you more exercises for your legs, feet and joints or you can find them on the (German-language) website www.my-medibook.de/uebungen-nach-diagnose.html (by Dr. Jens Wippert)

Simply scan the QR code for exercises.

Material composition

Bunions model • Hammer toe model:

90% cotton · 9% polyamide · 1% elastane

Pronation Control model:

45% polyester · 44% cotton · 9% polyamide · 2% elastane

Wash before use • wash inside out • do not use fabric softener

latexfrei
latex free
sans latex

Durability / Warranty

Store at room temperature and in a dry place. Avoid heat and moisture. If you observe all directions in this brochure the expected life-time will be 6 months maximum wearing time.

In case of item damage please contact your specialist retailer or dispensing chemist. Do not attempt to mend the taping socks yourself as this could affect the therapeutical effect.

The warranty period is 6 months starting from first use, as taping socks are subject to aging caused by material and wearing.

The warranty/guarantee also extends to the use of COMPRESSANA accessories.

Risks and side effects

Observe the instructions of the medical personnel. Taping socks can cause skin necroses, oedems and stress damage on peripheral nerves, in particular with inappropriate handling. If pain or negative changes (skin provoking) during ordered application should arise, immediately contact your doctor.

In case of sensory and/or wound healing disturbances (e.g. diabetic diseases) implicating a skin lesion, badly healing impairments can develop unnoticed. Patients of this group require special control as well as a specific instruction about particular risks and precaution possibilities.

If there are incompatibilities materials and substance of content please consult your doctor, pharmacist or specialist retailer before use.

The manufacturer is not responsible for damage and injuries, which developed due to inappropriate handling, misuse or wrong handling due to non-consideration of our references or use of accessories not originating from COMPRESSANA.

Disposal

No special disposal criteria exist. COMPRESSANA TAPE SOX can be disposed over the domestic refuse.

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COMPRESSANA GmbH